

STRADIS
HEALTHCARE

STRADIS
HEALTHCARE



Heart Cath Products

3025 Northwoods Parkway
Peachtree Corners, Georgia 30071
Telephone: (800) 886-7257
Fax: (770) 962-2391
www.stradishealthcare.com

Stradis Healthcare provides a wide variety of heart cath packs for your patient's needs. We serve more than 1500 surgical centers, private practices, and purchasing networks. Because we offer custom solutions, you can create the pack that's right for you.

The combinations are endless.

#686-070 CATH PACK

- 1 ea. Wrap, 45 x 45
- 1 ea. Drape, Bottom 76" x 76"
- 1 ea. Drape, Top 108" x 50"
- 2 ea. Drape, Side 44" x 76"
- 2 ea. Gown, Large, Non-Reinforced
- 1 ea. Bag, Header 16" x 22"
(10/cs)

#686-076 ANGIO PACK

- 1 ea. Cover Table, BTC 44" x 78"
- 1 ea. Drape, Femoral XL
- 1 ea. Bowl, Guidewire 2500CC
- 1 ea. Bowl, 32oz. 1000CC
- 1 ea. Bowl, 8OZ. 250CC, DENT
- 20 ea. Sponge All-Gauze, 4 x 4, 12-Ply
- 2 ea. Bag, Band Dome 50" x 35"
- 2 ea. Bag, Band 36" x 36" w/Tap
- 1 ea. Bag, Ziplock 9" x 12"
- 1 ea. Scalpel #11 Long Handle
- 1 ea. Skin Marker, Regular Tip
- 6 ea. Label, 1.5 X .5 WHT-CU

- 1 ea. Stopcock 3-way
- 2 ea. Syringe, 30CC L/L
- 4 ea. Syringe, 10CC L/L
- 1 ea. Scissors, Iris Curved
- 1 ea. Tipguard, Instrument
- 1 ea. Forcep, 5.5" Kelly
- 1 ea. Needle, 25G x 1.5"
- 1 ea. Needle, 22G x 1-1/2"
- 1 ea. Needle, 18G x 1.5
- 1 ea. Cover, Prove 6" x 48"
- 2 ea. Conductivity Gel, 20ML
- 1 ea. Chloraprep, 2% 10.5ML
- 2 ea. Rubber Bands
- 1 ea. Needle Counter 10CT
- 1 ea. Fluid, Management System
- 1 ea. Drape, C-Arm 41 x 74
- 6 ea. Towel, Blue O.R.
- 1 ea. Drape, Set Up. Clear
- 1 ea. Gown, Large, Non-Reinforced
- 1 ea. Gown, XL, Non-Reinforced
- 2 ea. ABS Towel 15" x 26"
- 1 ea. Bag, Header 19" x 30"
(3/CS)



CALL NOW FOR YOUR FREE CUSTOM SAMPLE.
CALL 1.800.886.7257 | WWW.STRADISHEALTHCARE.COM